

Analytical Balance Internal Calibration



- Attractive Modern Design .
- Constructed with Cast-Aluminium external body.
- The panel with five operation button, sliding glass windscreen provide big space, full vision and easy operation .
- Use high precision magnet sensor to achieve fast and stable reading.

Specification of Analytical Balance Internal Calibration

Model	ST-FA1004N	ST-FA1104N	ST-FA1204N	ST-FA2004N	ST-FA2104N	ST-FA2204N
Capacity (g)	100g	110g	120g	200g	210g	220g
Resolution (g)	0.0001g					
Min weighing (g)	0.0004g					
Stable Time	≤ 3 S					
Display	LCD (white back light with black font)					
Pan Size	Ø80mm					
Operation temp	5-25°C					
Repeat ability	±0.0002g					
Liner	±0.0003g					
Draft shield size	240 x 190 x 265 mm					
Cal.weight	Internal Calibration					
Option	RS232C & RS485, Printer					
Packing Size	Inside: 415 x 285 x 430mm			Outside: 475 x 315 x 485mm		
G.weight	10Kg					