



THREE LEVEL MICROTUBE RACK FOR 0.2 - 0.5 AND 1.5 ml

Made from polypropylene and can be linked at either end to form longer racks. Numbered positions. Lower level 24 micro tubes of 0.2 ml. Middle level 14 micro tubes of 0.5 ml. Upper level 12 micro tubes of 1.5 ml. Standard colour yellow. Part No. **1001542**



MULTI SUPPORT TUBE RACK TUBE SIZE Ø 10/13 AND 14/17 mm

Suitable for drying and storing of tubes. Autoclavable. Suitable for use in a water bath. Standard colour blue. Part No. **1001543** Group of 80 spikes for tubes of 10-13 mm Ø. Part No. **1001544** Group of 50 spikes for tubes of 14-17 mm Ø.



REVERSIBLE TUBE RACK FOR MICROTUBES 0.2/0.5 AND 1.5/2 ml

Made from polypropylene with removable lid that fits on any side. Alphabetic identification positions. Face A for 168 micro tubes of 0.2 ml and 12 x 1.5 ml. Face B for 110 micro tubes of 0.5 ml and 12 x 1.5 ml. Standard colour yellow. Part No. **1001545**



REVERSIBLE TUBE RACKS FOR 0.5 - 1.5 AND 2 ml TUBES

Face A for 96 micro tubes of 0.5 ml. Face B for 96 micro tubes of 1.5 to 2 ml. Made from polypropylene with removable lid. Alphabetic tube position identification. Height 26 mm, with lid 50 mm. Width 115 mm. Length 250 mm. Standard colour yellow. Part No. **1001524**



BOX TO STORE 100 CRYO VIALS AND 1.5 TO 2 ml MICROTUBES

Made from polyethylene, with hinged, lid and pressure latch close. Alphabetical grid tube location within the lid. Height with lid 55 mm. Width 140 mm. Length 140 mm. Standard neutral colour. Part No. **1001525**



BOX FOR 50 CRYO VIALS AND 1.5 TO 2 ml MICROTUBES

Made from polyethylene, with hinged, lid and pressure latch close. Alphabetical grid tube location within the lid. Height with lid 55 mm. Width 80 mm. Length 140 mm. Standard neutral colour. Part No. **1001526**



STORAGE BOX FOR 96 TUBES OF 0.2 ml

Made from polypropylene, with hinged lid and pressure latch close. Alphanumeric identification positions. Height with lid 30 mm. Width 90 mm. Length 130 mm. Standard colour yellow. Part No. **1001527**



BOX FOR STORING 12 X 10 mm CUVETTES

Made from Polypropylene, with hinged lid and pressure latch close. Removable base. Cuvette heights up to 60 mm high. Height with lid 60 mm. Width 115 mm. Length 115 mm. Standard colour neutral. Part No. **1001529**



FREEZER BOXES FOR VIAL STORAGE

Made from polycarbonate. Suitable for freezer or liquid nitrogen storage for temperatures between -196 °C to +121 °C. Alphanumeric identification, with ventilation holes and drain. Pressure lid. Part No. **1001713** Box of 25 vials. Height with lid 55mm. Width 75 mm. Length 75 mm. Part No. **1001714** Box of 100 vials, Lid with hinge. Height with lid 55mm, Width 130mm. Length 130 mm.