



MULTI PURPOSE KNIFE

With a black bakelite® handle. Designed for comfort and safety
 Made of tempered stainless steel.
 Total length: 230 mm.
 Blade length: 95 mm.
 Part No. **1000862**



VIBRATING SPATULA AND BALANCE WEIGHT FORCEPS

Vibrating spatula.

With manual vibrating device. Made of AISI 304 stainless steel.

Part No. **1000413** Curve end. Total length 176 mm. Width 12 mm.

Part No. **1000412** Grooved end. Total length 176 mm. Width 12 mm.

Balance weight forceps.

Made of AISI 304 stainless, electrobright finish. Total length 93 mm.

Part No. **1000217**

Forceps for balance weights.

Made of AISI 304 stainless, electrobright finish. Total length 100 mm.

Part No. **1000216**



SPATULA POCKET POUCH

Part No. **1001530**

Contents:

1 Flexible spatula. Part No. 1000856 . 1 Semi micro spatula. Part No. 1278140

1 Micro spatula flat spoon. Part No. 1277150. 1 Dissection scissors. Part No. 1000807

1 Forceps blunt end. Part No. 1000580



SODIUM SPOON

With brass basket. Ø 25 mm. Total length: 190 mm.

Part No. **1000805**



FUSION SPOONS

Made of AISI 304 stainless, wooden handle. Total length: 350 mm.

Flat base. Volume 90 ml.

Part No. **600093**

Round base. Volume 75 ml.

Part No. **600094**



POWDER SCOOPS

Length: 100 mm. Width: 26 mm.

Part No. **6000278** Stainless steel AISI 304.

Part No. **6001278** Pure nickel.



SCOOPS WITH HANDLE

Made of stainless steel AISI 304.

Part No. **1001592** Medium Total length 320 mm.

Part No. **1001593** Big Total length 360 mm.



SCOOPS WITH HANDLE

Made of stainless steel AISI 304.

Part No.	Total length	Usable length	Width	Capacity
1000859	180 mm.	100 mm.	50 mm.	120 ml.
1000860	260 mm.	150 mm.	75 mm.	450 ml.
1000861	310 mm.	200 mm.	100 mm.	1500 ml.
1001354	330 mm.	210 mm.	150 mm.	2000 ml.