

Sastec Mobile Clean Sampling Cabinet



Features:

- Constructed in accordance to GMP Standard, with cleanliness to Class 1000.
- Three sides with polyglass material and entrance with flexible strips on lockable castors.
- Top box consists of HEPA Filter with centrifuged blower to deliver 99.995% of clean air against 0.3um of particles.
- With fluorescent Light, control built at the entry support stand including digital airflow adjustable from 0.36 to 0.45 m/s

Specification of SASTEC Mobile Clean Sampling Cabinet	
Model	ST-MCC3
External Dimension	740 x 740 x 1950 mm
Internal Dimension	600 x 650 x 1550 mm
Material	All Stainless Steel
Side Panels	Acrylic Polyglass
Voltage	220 - 230V 50Hz
Gross Weight	110Kg
Packing Size	850 x 850 x 2050 mm