



ROTATORY EVAPORATOR "RS 3000-V"

FOR 1000, 2000 AND 3000 ML. FLASKS.
MULTIFUNCTIONAL: VACUUM PUMP AND BATH BUILT-IN

Quality and safety. Rotary speed and bath temperature regulation

FEATURES

Compact rotatory evaporator composed of stainless steel AISI 304 case and with an extra-soft manual elevator system and bath made of dur-al alloy chemically treated.

Rotatory speed: from 20 to 120 rpm.

Bath temperature: up to 100 °C.

External probe temperature: up to 100 °C.

Vacuum breaker with air inlet filter.

CONTROL PANEL



1. Starting main switch ON/OFF.
 2. Digital reading indicator of external probe temperature.
 3. Bath temperature digital regulator.
 4. Rotator pump switch.
 5. Heater bath switch.
 6. Water recirculation pump switch.
 7. Regulator knob of rotating speed.
- Vacuometer built-in the case.



MODEL

	Part no	Speed Regulation t.p.m.	Height / Width / Depth (extérieur) cm			Power W	Weight Kg
RS 3000-V	6003000	0 a 120	64	74	42	1050	38

It is completely supplied with a tap for solution, a PT100 probe, a coolant, two 3000 ml. flasks, one 2000 ml flask, two 1000 ml. flasks and its accessories.

GLASS AND ITS ACCESSORIES DIAGRAM SPARE PARTS

