

ST-L500 Tabletop Low Speed Centrifuge



- ST-L500 Tabletop Low Speed Centrifuge is designed for the clinical research.
- Stainless rotors with variety of capacity for convenient operation.

Rotor Options



Order No.: 808350
Max. Speed: 5000rpm
Max. RCF: 4030xg
Rotor Capacity: 24 x 10ml
Adaptor: 808308 (1.5ml)

Features

- Microprocessor control and DC brushless motor.
- Touch panel and LCD display.
- RCF value can be calculated automatically.
- Special damping structure to reduce vibration.
- Electrical door interlock, with which centrifuge cannot work if the door is open and the door cannot be opened when it is working.
- Variety of brackets for users to convenient



Order No.: 855201
Max. Speed: 4200rpm
Max. RCF: 2760xg

Bracket



Order No.: 855210
Capacity: 4 x 50ml
Adaptor: 854708 (15ml falcon)



Order No.: 855220
Capacity: 4x4x15ml
Adaptor: 854806 (1.5ml falcon),
854808 (5ml falcon)



Order No.: 855230
Capacity: 4x6x5ml
Adaptor: 854809(7ml), 854807(1.5ml)
854813(0.5ml)

Specification of Tabletop Low Speed Centrifuge

Model	ST-L500
Order No.	209000
Max. Speed	5000rpm
Max. RCF	4030xg
Max.Capacity	16 x 15ml
Speed Accuracy	±30rpm
Time Setting Range	1min to 99min
Noise	≤65dB(A)
Power Supply	AC220V±22V 50/60Hz 5A
Wattage	250W
Dimension (W x D x H)	380 x 500 x 280mm
Package Dimension (W x D x H)	660 x 465 x 350mm
Net Weight	19kg