

Water Bath Economic Type



- Interior material is made by one-step forming technology to keep the entire liner free from any joints or welds, so as to eliminate the problems of rust, corrosion and water leakage.
- Microprocessor PID temperature controller with temperature control protection, digital display and timing function ensures a precise reliable temperature control.
- Option for stainless steel cover lid or reduction ring set.

Specification of Water Bath

Model	ST-WBB21	ST-WBB22	ST-WBB24	ST-WBB21
Internal Dimension (W*D*H) mm	160 x 160 x 130	300 x 150 x 140	300 x 300 x 140	450 x 300 x 130
Capacity	3 Litre	6 Litre	12 Litre	17 Litre
Heater	300W	600W	800W	1500W
Note	One holes	Two holes	Double Columns Four Holes	Double Columns Six Holes
Material	Internal	Stainless Steel (one step forming body)		
	External	Power Coated Steel		
Temperature	Range	Ambiet+5 ~ 99°C		
	Accuracy	±0.1 °C		
	Uniformity	0.5°C		
	Sensor	PT100		
	Controller	Digital Control LED Display		
Timer	0~ 99 min			
Power Supply	AC 110V~ AC 127V; AC 220V~ AC 240V; 50Hz ~ 60Hz;			