



Hotplate-Rectangular "Combiplac"

PLATE TEMPERAURE ADJUSTABLE UP TO 400 °C.

FEATURES

Electronic energy temperature control.
The non warping steel hotplate has heating elements embedded throughout the hotplate surface maintaining the heat distribution. 
The main case is made of AISI 304 stainless steel.

CONTROL PANEL

Main power switch with "ON" indicator lamp.
Electronic energy control from 0 to 100 %.
An indicator lamp shows when the unit is hot.
Connector at the back for the external digital thermometer, "Sensoterm" or "Electemp".
See Regulators and Controllers (pages 286 and 287).

MODELS

Part No.	Width / Depth (plate) cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
3000156	20 40	12 22 55	2300	9
3000718	30 30	13 33 46	2800	13



ACCESSORIES

Sand bath adapters. Made entirely of AISI 310 stainless steel, removable for easy cleaning.

Part No.	Hotplate Part No.	Height/Width/Depth (Sand tray dims) cm
6000719	3000156	5.5 20 40
6000720	3000718	5.5 31 31



Hotplate-Circular "Combiplac"

PLATE TEMPERATURE ADJUSTABLE UP TO 400 °C.

FEATURES

Electronic energy controller and embedded heating element throughout the non warping steel hotplate surface. 
A retort bar can be located in a hole at the back of the unit.
The main case is made of AISI 304 stainless steel.

CONTROL PANEL

Main power switch with "ON" indicator lamp.
Temperature controlled by a 0 - 100 % energy controller.
A lamp indicates when the hot plate is "ON".
A Thermometer "Sensoterm" or "Electemp" that can be added and are listed as accessories.
See section Regulators and Controllers (pages 286 and 287).



MODELS

Part No.	Plate Ø cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
3002301	12	10 16 27	550	2
3002302	19	13 21 30	950	3

ACCESSORIES

Suitable for the 12 cm Ø hotplates only Part No.3002301:

Bath of AISI 304 stainless steel of 140 mm Ø x 75mm high, with concentric reduction ring lid, Suitable for water, oil or sand. Part No. **7002510**

Round bottom flask adapter made of aluminium and suitable for Kjeldahl flasks of 100 ml and 50 to 100 ml round bottom flasks. Part No. **1002290**

Round bottom flask adapter made of aluminium and suitable for Kjeldahl flasks of 300 to 500 ml and round bottom flasks of 250 ml. Part No. **1002280**

Support rod of 12mm Ø and 500 mm high, fits on to the main body of the hotplate. Suitable for the circular Combiplacs 3002301 and 3002302. Part No. **1000725**



Electronic contact thermometers with digital readout "Sensoterm". With Pt100 sensor with stainless steel AISI 304 sheath. Part No. **1001009**

Features. See Regulators and Controllers (pages 286 and 287).

Electronic temperature controller "Electemp"

Part No.	Height / Width / Depth (exterior) cm	Max power (at 230 V) W	Weight Kg
3000887	10.5 13.5 21	3100	1.250

Features. See Regulators and Controllers (pages 286 and 287).

