



Stainless Steel Rack-Double Sided with Locking Wheels



Specifications:

- Made in AISI 304 stainless steel with 4B finish
- Rack to accept high and low baskets with measurement of 585 x 395 x 195 mm and 585 x 395 x 100 mm
- With 4 locking castors
- Dimensions: 1800L x 850D x 1780H mm
(Other dimensions are available upon request)



Loading & Unloading Trolley For Steriliser With Rack

Specifications:

- Made in AISI 304 stainless steel with 4B finish
- It has loading bracket for upper rack to slide into the chamber
- With adjustable height to allow trolley to be used with most Steam Sterilizers.
- Upper rack has stainless steel roller
- Dimensions of trolley: 1600L x 600W x 700-800H mm
- Dimensions of racks: 1600L x 600W x 600H mm
(Other dimensions are available upon request)



Stainless Steel Baskets



Specifications:

- Made in AISI 304 stainless steel, electro-polished surface.
- Material with 3mm body and 5mm surround.
- Dimensions: 585 x 395 x 100 mm
- 585 x 395 x 195 mm
(Other dimensions are available upon request)