

ST-JA Electronic Balance

- Cast-aluminium outer covering. High weighing value and precise balance. Easy to operate, stable rapidly with taring and auto zero tracking function. Serial interface RS232/USB



Specification of Electronic Balance

Model	ST-JA3000AC	ST-JA4000AC	ST-JA5000AC
Capacity (g)	3100g	4100g	5100g
Sensor		Load Cell	
Resolution (g)		0.01g	
Min Weighing (g)		0.04g	
Stable Time		≈ 3S	
Display		LED/ LCD	
Pan Size		180 x 180mm	
Operation Temp		5-25°C	
Repeat Ability		±0.03g	
Cal. Weight		2000g	Option
Option		RS232C & RS485 , Printer, USB	
G. Weight		5.0Kg	

** The calibration weight is optional