

ST-XY-B Electronic Balance

- With digital LED display ST-XY-B series is widely used in industrial applications, laboratory and university under weighing functions to test the density is included.



Specification of ST-XY Electronic Balance

Model	ST-XY600B	ST-XY1000B	ST-XY2000B	ST-XY3000BF	ST-XY5000BF
Capacity (g)	610g	1100g	2100g	3100g	510g
Resolution (g)	0.1g				
Min. Weighing (g)	0.4g				
Stable Time	≤2S				
Pan Size	Ø130mm			160 x 160mm	
Operation Temp	5-25°C				
Repeat Ability	0.2g				
Liner	0.3g				
Cal. Weight	200g	500g	1000g	2000g	
Option	RS232C & RS485, Printer				
Packing Size	Inside: 310 x 230 x 130mm Outside: 660 x 470 x 320mm				
G. Weight	2.0/ 22Kg				

** The calibration weight is optional