

Upright Ultra Low Freezer



ST-UR40-270

ST-UR40-450

- Temperature Range -10°C to -40°C.
- Microprocessor temperature controller with LCD display.
- Display accuracy $\pm 0.1^\circ\text{C}$ with Pt 1000 sensor.
- Safety door locking system.
- External is made of quality steel, internal is ABS material.
- Audio/Visual alarm system for high/low temperature.
- 110mm CFC-Free insulation.
- Supplied with 7 drawers.
- CFC Free refrigerant R507
- Certification: ISO 9001, ISO 14001, ISO 13485
- Power supply: 220V 50Hz (110V 60Hz is available)
- Warranty period: 2 years on all parts.

Optional:

- 7 days Chart Recorder
- 5 years extended warranty

Specification of Upright Ultra Low Freezer

Model	ST-UR40-270	ST-UR40-450
Capacity (L)	270L	450L
External Dimension (W x D x H) mm	700 x 600 x 1750	660 x 780 x 1850
Internal Dimension (W x D x H) mm	440 x 400 x 1490	
Input Power (W)	370W	
Weight (kg)	70Kg	180Kg
	Supply with 7 drawers	12 drawers audible/ visual for: High/ Low temperature sensor failure and system failure