



Colorimeter "WSD-3"

THE WSD-3 IS A HIGH PRECISION SPECTRAL COLORIMETER FOR MEASURING THE COMPOSITION OF THE COLOUR IN CHROMATIC COMPONENTS. THE OPTICAL SYSTEM WORKS IN A REFLECTION MODE TO EVALUATE THE CHROMATIC COMPONENTS OF ANY SURFACE.



APPLICATIONS

Chromatic, surface composition analyser.
Suitable for research and quality control.

FEATURES

- Illumination/measuring conditions:o/d
- Measuring system:reflection
- Spectral range400 to 700 nm
- Illumination standard:D 65
- Measuring time:10 sec.
- Light source Halogen lamp:6V 12W
- Head entrance diameter:22 / 6 mm
- Field of vision10°
- Reading precision:X,Y ± 0.0001, other scales ± 0.01
- Measuring error:0.1
- Stability {Zero set ≤ 0.1
.....Measuring error ≤ 0.1
.....VPower variation ≤ 0.1
- Measuring precision:.....0.1
- Preheat time:10 min.
- Ambient conditions:0 - 40 °C
- Power requirement:.....220V 50Hz 10%



MODEL

Part No.	Height / Width / Depth cm	Power W	Weight Kg
5120190	29 30 32	23	4.5

SPARE

Halogen lamp 6V/12W.
Part No. **5313040**



Microtiter Strip Reader "Reader M-2000"

STORAGE MEMORY FOR UP TO 59 METHODES INCLUDING: FACTORS, STANDARDS OR CALIBRATION CURVES WITH UP TO 10 STANDARDS. ADDITIONAL CALCULATION FEATURES.

APPLICATIONS

Clinical diagnostics, E.I.A., Veterinary, quality control.

FEATURES

Spectral range: 400 to 700 nm, using interference filters of 10 nm band pass.
Supplied complete with a filter of 450 nm (other wavelengths available consult for options.)
Measuring range: -0.2 to 3.5 O.D.
Photometric accuracy: >1%.
Photometric precision: ±1%.
Light source: Long life LED.

Detector: solid state.
59 method memories.
Curve point editing feature.
Software: ABS. Concentration and polynomial curve fitting with up to 10 standards.
RS-232 output to connect a printer or to a PC.
Interface RS-232.



MODEL

Part No.	Height / Width / Depth cm	Power W	Weight Kg
4120013	11 18 28	7	2

ACCESSORY

Ink printer (not thermal). Paper 2 1/4" (56 mm) wide roll. Interface RS232. Includes power transformer and interface cables. (5V-3A)



Model

Part No.	Height / Width / Depth cm	Weight Kg
4120113	4 16 10	1

SPARE

Lamp of 4 mm Ø 5 V.
Part No. **4313030**