

Liquid Nitrogen Tank Accessories

Accessories For Liquid Nitrogen Tank	
Model	Description
ST-C430	Cart for liquid nitrogen tank size 30-50L
ST-C520	Cart for liquid nitrogen tank size 100L and above
ST-CAN4	Square canister with 4 shelves and plastic box
ST-CAN5	Square canister with 5 shelves and plastic box
ST-CAN19	Round canister 19mm
ST-CAN26	Round canister 26mm
ST-CAN38	Round canister 38mm
ST-CAN63	Round canister 63mm
ST-CAN94	Round canister 94mm

Differential Blood Counters ST-DBC-9

- Used in hematology to classify and count white blood cells.
- Outer shell resist corrosion.
- At one hundred increment of the total record, it will ring.
- Reset knob will clear all figures.
- Reading 0~999.
- Dimension (L x W x H): 320 x 80 x 50mm

