


Rotary Evaporator, LCD Digital (ST-RE100-PRO)



Features:

- Easy to read large digital LCD screen displays heating temperature, rotation speed and timing.
- Speed range from 20 to 280 rpm
- Large 5L water-oil heating bath with heating temperature range of RT to 180°C.
- Heating bath with precise temperature control and adjustable safety circuit.
- Designed a detachable operating panel allow remote control.
- Patented condenser (cooling surface 1,500 cm²) with excellent cooling effect.
- Evaporating flask up to 3 litre.
- Vertical and diagonal glassware set is available.
- Motorized lift with quick-action, automatic release evaporating flask to top position in case of power failure.
- Adjustable final position recognition to protect operator and sample against breakage.
- Evaporating flask with an ejector, convenient to remove.
- Available with timer function to precise control processing.
- Chemical resistant double PTFE system and patented pressure spring provide excellent sealing.
- PC can control and document all parameters via RS232 interface (USB).