

LCD Display



Specification of Rotary Evaporator LCD Digital

Model	ST-RE100-PRO
Voltage	100-120/200-240
Frequency	50/60 Hz
Power	1400W
Motor Principle	DC brushless motor
Speed Range	20-280 rpm
Speed display	LCD
Reversible Direction of Rotation	Yes
Heating Temperature Range	180°C
Heat Control Accuracy	±1°C
Heat Output	1300W
Temperature Display	LCD
Lift	Motor
Stroke	150 mm
Timer	Yes
Timer Display	LCD
Time Setting Range	1-999 min
Overall Dimension (D x W x H) mm	465 x 457 x 583
Weight	15 kg
Permissible Ambient Temperature	5-40°C
Permissible relative Humidity	80%
Protection Class acc. to DIN EN60529	IP20
RS-232 Interface	Yes