


Magnetic stirrer



ST-MSD

ST-MSA

- Compact design.
- Stainless steel surface.
- Infinitely variable speed motor.
- Motor with electronic speed control, constant speed even during changes in load.
- High Magnetic adhesion.
- Enclosed assembly (IP42) guarantees long service life.
- Easy to read backlit LCD display for digital model.
- CE Marking.

Specification of Magnetic Stirrer

Model	ST-MSA (Analog)	ST-MSD (Digital)
Power		20W
Maximum Stirring Capacity (HzO)		20 Litres
Maximum Stirring Bar Length		80mm
Motor Type		Brushless Motor
Motor Input Rating		19.2W
Motor Output Rating		10W
Speed Range		0 ~ 1500 rpm
Material		Stainless Steel Plate, Ø 135mm
Dimension (W x D x H) mm		200 x 130 x 65mm
Weight		2.4kg
Permissible Ambient Temperature		5 ~ 40°C
Permissible Related Humidity		80%
Protection Class		IP42 (EN60529)
Voltage		85 ~ 265V, 50/60Hz