

## Drying Oven, Stainless Steel



- Both internal and external made of non-rusting Stainless Steel.
- Adjustable Shelves.
- Auto Controller of fan speed.
- Pre-diagnostic function.
- Over temperature alarm.

### Optional

- Independent audible and visible temperature alarm.
- Programmable Controller with 15 segments, timer 1 to 99 hour for each segment and fan speed 0 to 100% adjustable.
- RS 485 for connection to computer and printer.

### Specification of Drying Oven

Model	ST-CR31	ST-CR51	ST-CR91	ST-CR32	ST-CR52	ST-CR92
Temperature Range	10°C above ambient ~200°C					
Display Resolution	0.1°C					
Temperature Stability	±1°C					
Temperature Uniformity	2.5%					
Power Consumption	200W	400W	600W	200W	400W	600W
Electrical Requirement	220V 50Hz					
External Dimension (mm)	600 x 500 x 520	650 x 600 x 600	700 x 660 x 740	600 x 500 x 520	650 x 600 x 600	700 x 660 x 740
Internal Dimension (mm)	300 x 310 x 320	350 x 360 x 400	400 x 420 x 536	260 x 280 x 330	310 x 330 x 480	400 x 430 x 520
Capacity	24L	50L	90L	24L	50L	90L
Shelves (Pcs)	1	1	2	1	1	2
Weight (Kg)	23 Kg	30 Kg	41 Kg	23 Kg	30 Kg	41 Kg
Air Circulation	Natural	Natural	Natural	Fan	Fan	Fan
	2310.00 / 2772.00	3110.00 / 3732.00	4370.00 / 5244.00	2900.00 / 3480.00	3580.00 / 4300.00	4980.00 / 5980.00