



Chromatography tanks

APPLICATIONS

For quantitative analysis. Useful for storing and protecting the plates by means of the hermetic upper lid.

FEATURES

Made of heavy thick glass, with a flat bottom, ground rims with rounded edges. Solid glass rectangular lid, smooth base with ground edges to make an air tight fit tank. Rectangular glass lid with handle, smooth base and ground surface for hermetically sealed closure with the tank.

MODELS	Part No.	Dimensions mm	N° of plates	With Lid
MINICROM	1001162	100 x 100	1	Yes
CROMACROM	1001163	200 x 200	1	Yes
MULTICROM	1001164	200 x 200	5	Yes

The model MULTICROM has 5 plate slots to maintain the plate upright.



Stackable Anatomy - Histology sample storage chest

FEATURES

Manufactured in steel epoxy coated. Resistant to corrosive atmospheres.

STACKABLE MODULAR SYSTEM

Cabinet for storing slides.

Complete with 14 drawers for slides of 26 x 76 mm.

Drawer capacity: 350 slides.

Part No. **1052014**

Cabinet for storing transparencies.

Complete with 7 drawers for transparencies of 50 x 60 mm.

Drawer capacity: 230 transparencies.

Part No. **1005207**

Cabinet for storing paraffin wax mounted samples.

Complete with 10 drawers.

Part No. **1052010**

MODELS

Part No.	Height / Width / Depth (usable) cm	Height / Width / Depth (exterior) cm	Weight Kg
1052014	9 2.6 45	14 48 48	9.5
1005207	7.2 5.8 47	12 12 48	6.5
1052010	2.6 43.5 47	34.5 48 48	18.5

Base and lid units

Mounted on the top or as a support base. Fits the individual modules to form distinct user definable cabinets.

Part No.		Height / Width / Depth (external) cm	Weight Kg
1005208	Base	5,5 48 48	2,5
1005209	Lid	5,5 48 48	2,5

