

Portable pH Meter



- ST-PH4- Portable pH Meter with case and electrode

- ST-PH260- Portable pH Meter with case software and electrode

Specification of SASTEC Portable pH Meter

Model	ST-PH260	ST-PH4
Parameters	pH/ mV (ORP)°C	pH/ mV (ORP)°C
Range	0.00 ~ 14.00	
pH		
mV	± 1800	± 1400
Temperature	-5.0 ~ 105.0°C	-
Resolution	0.01	
pH		
mV	1	
Temperature	0.1°C	-
Accuracy	± 0.01	
pH		± 0.03
mV	± 0.1%	± 0.2%
Temperature	± 0.3°C	-
Temperature Compensation	Manual/ Auto: -5.0 ~ 105.0°C	Manual: 0.0 ~ 60.0°C
Calibration	May select 5 standard buffer solution for one and Two-point calibration	Two-point calibration
Data Storage	200 data with time-and-date	-
View, Delete and Print	Yes	-
Display	Backlit LCD	
Under-Voltage Indication	Yes	-
Output	RS-232	-
Software	1.1DC	-
Housing	IP65	-
Power	Four AA Batteries	Two AA Batteris
Case Outside Serial No.	WXS-A006-1	WXS-A004-1
Dimension (mm)	210 x 86 x 50	170 x 75 x 30
Weight (Kg)	0.5 Kg	