



## Water bath with inner stirring “Agibat-20” with digital control and display of temperature and time

ADJUSTABLE TEMPERATURE RANGE FROM AMBIENT +5°C UP TO 80°C. STABILITY  $\pm 0.1^\circ\text{C}$ .  
SET TEMPERATURE ACCURACY  $\pm 1^\circ\text{C}$ , WORKING TEMPERATURE PRECISION  $0.1^\circ\text{C}$ .

NEW

*P.Selecta original model. More space. More accuracy. Economical.*



*Stirring detail by internal circulation pump.*

### FEATURES

Double skinned tank, inner case made of AISI 310 stainless steel.  
External case of AISI 304 steel.  
Submerged heating elements in INCOLOY special alloy.  
It comes complete with drain tap.

### STIRRING SYSTEM

These baths have an aspiration and outflow internal circulation pump which allow a constant flow stirring.  
Flow: 8 litres per minute.

### CONTROL PANEL

See page 81.  
Three positions switch: only heater, heater + stirring and start/stop.

### MODEL WITH MAKROLON® FOLDING LID

Due to its transparency, you can visualize the bath's inside.  
It improves temperature stability, and restricts evaporation and pollution.

### MODELS

Part No.	Lid Makrolon®	Maximum Temp °C	Litres Capacity	Height/Width/Depth (usable) cm	Height/Width/Depth (exterior) cm	Power W	Weight Kg
6032012	Yes	80	20	15 48 30	20 58 42	2010	11
6032013	No	80	20	15 48 30	40 58 42	2010	12

All baths with heating shelf and drainage. Faucet included.

### ACCESSORIES.

See page 83.

### SAFETY:

DIN 12877.2 STANDARD. THE UNIT WILL CUT OUT IF THERE IS INSUFFICIENT LIQUID OR IF THE BATH TEMPERATURE EXCEEDS THE SET TEMPERATURE. A MANUAL RESET ADJUSTABLE TEMPERATURE SAFETY THERMOSTAT IS PLACED AT THE BACK OF THE UNIT.

