

## Laboratory Refrigerator



**ST-RG-260**



**ST-RG-950 (2 Door)**



**ST-RG-1500 (3 Door)**

- Temperature Range +2° to 10°C, can be set freely.
  - Digital temperature controller.
  - Pt 1000 sensor with 0.1° accuracy on the display.
  - 2 - Layer glass door
  - Forced air cooling with large ventilator.
  - Audio/visual alarm for high/low temperature, door ajar, and system failure.
  - Interior with fluorescent lighting.
  - 5 wire shelves.
  - Exterior made of high quality steel, interior is ABS material.
  - CFC-Free R134a refrigerant.
  - Mounted on castors.
  - Power supply: 220V 50Hz (110V 60Hz is available)
  - Certificate: ISO 9001, ISO 14001, ISO 13485
  - Warranty period: 2 years on all parts.
- **Optional:** i. Printer                      ii. 5 years extended warranty

### Specification of Laboratory Refrigerator

Model	ST-RG-260	ST-RG-300	ST-RG-950	ST-RG-1500
Capacity (L)	260	300	950	1500
External Dimension (D x W x H) mm	590 x 603 x 1556	590 x 603 x 1728	780 x 1200 x 1894	800 x 1600 x 1900
Internal Dimension (D x W x H) mm	445 x 525 x 1070	445 x 530 x 1265	600 x 1100 x 1400	650 x 1500 x 1600
Compressor HP	1/5	1/3	3/8	1
Weight (kg)	105	110	180	220