

## High Speed Mini Centrifuge



- Powerful centrifuge force 15100 g, broad application range
- User friendly large LCD display provides all visual information
- Brushless motor drives quickly and effortlessly accelerates the rotor to set speed
- Gravity centre strutual design, high strength metal casing, make sure stabilization and security
- Compact design, save valuable space
- High speed, low noise
- Fully autoclaveble high strength aluminium alloy rotor
- High performance price ratio and cost effective operation

Specification of High Speed Mini Centrifuge	
Model	ST-D2012
Max. Speed	15,000 rpm (500-15,000 rpm), increment: 100 rpm
Max. RCF	15,100 xg, increment: 100 rpm x g
Rotor	2ml x 12
Run Time	30 sec ~ 99min
Driving Motor	Brushless DC Motor
Safety Devices	Door Interlock, Over Speed detector, Over Temperature Detector, Automatic Internal Diagnosis
Power	100V-240V 50/60Hz, Single Phase
Dimension (L x D x H)	255 x 245 x 140mm
Weight	6 kg
Additional Features	Speed/Acceleration Switch, Short time run function, sound alert function