

Electronic Balance FA/ JA Series

- Modular Electromagnetic force sensor.
- White light LCD Display.
- High sensitivity touch buttons.
- Under hook weighing.
- The glass door can be locked during transportation.
- RS232 interface module, USB interface module, built-in clock.
- Independent "Tare" and "Zero" button, two-way communication
- Four level Shockproof
- Dynamic temperature compensation
- Weighing speed adjustable, showing way adjustable



ST-FA224

Specification of Electronic Balance FA Series

Model	ST-FA124	ST-FA224	ST-FA1004	ST-FA1104	ST-FA1604	ST-FA2004	ST-FA2104	ST-FA2104	
Weighing Capacity (g)	120g	220g	100g	110g	160g	200g	210g	60/ 210g	
Readability (mg)	0.1mg							0.1/ 1mg	
Repeatability (≤mg)	±0.1 ≤mg							±0.1/ ±1 ≤mg	
Linearity (≤mg)	±0.2 ≤mg							±0.2/ ±2 ≤mg	
Pan Size (mm)	Ø 90mm								
Gross Weight (Kg)	7.5Kg								
Net Weight (Kg)	5.5Kg								
Dimension (mm)	365 x 223 x 338mm								
Packing Measurement (mm)	500 x 310 x 450mm								

Specification of Electronic Balance JA Series

Model	ST-JA1003	ST-JA1203	ST-JA2003	ST-JA3003	ST-JA4103	ST-JA5003
Weighing Capacity (g)	100g	120g	200g	300g	410g	500g
Readability (mg)	0.001g					
Repeatability (≤mg)	±0.001 ≤g					
Linearity (≤mg)	±0.002 ≤g					
Pan Size (mm)	Ø 90mm					
Gross Weight (Kg)	7.5Kg					
Net Weight (Kg)	5.5Kg					
Dimension (mm)	365 x 223 x 338mm					
Packing Measurement (mm)	500 x 310 x 450mm					