

Upright Fridge-Freezer

- Temperature Range- Fridge : +2°C to 8°C
Freezer: -10°C to -40°C
- Microprocessor temperature controller, LCD display
- Audio/Visual alarm for High/Low temperature
- 90mm CFC-free insulation
- 2 units of compressor (Danfoss) and EBM motor fan.
- Double sealing with lock
- Painted metal external and ABS material internal.
- Three shelves for upper refrigerator compartment and 3 drawers for lower freezer.
- Mounted on castors.
- CFC-free refrigerant R507 and R134a.
- Certification: ISO 9001, ISO 14001, ISO 13485
- Power supply: 220V 50Hz (110V 60Hz is available)
- Warranty period: 2 years on all parts.



ST-FF40-208

Specification of Up Right Fridge-Freezer

Model	ST-FF40-208	
Type	Refrigerator	Freezer
Capacity (L)	154	99
External Dimension (W x D x H) mm	623 x 604 x 1775	
Internal Dimension (W x D x H) mm	397 x 368 x 781	397 x 368 x 414
Input Power (W)	400	
Weight (kg)	75	