

## Accessories


### Glassware

Compact vertical condenser for all standard distillations.

Condenser type: Vertical

Cooling Surface: 1500cm<sup>2</sup>

Catalogue No: 18900165



Evaporating Flask NS 29/32	
Catalogue No.	Volume
18300116	50ml
18300117	100ml
18300118	250ml
18300119	500ml
18300120	1000ml
18300121	2000ml
18300122	3000ml



Receiving Flask KS 35/20	
Catalogue No.	Volume
18300107	100ml
18300108	250ml
18300109	500ml
18300110	1000ml
18300111	2000ml
18300112	3000ml




Spiders with 5 flask NS 29/32	
Catalogue No.	Volume
18300102	50ml
18300103	100ml



Connector NS 29/32	
Catalogue No.	Size
18300094	NS 29 / NS 24
18300095	NS 29 / NS 19
18300096	NS 29 / NS 14



Foam Brake, 250ml NS 29/32  
Catalogue No: 18300093



Vapor Tube NS 29/32  
Catalogue No: 18300097



Seal  
Catalogue No: 18300098

### Vacuum Pump Catalogue No: 18200743



Voltage (V)	110/220	
Frequency (Hz)	50/60	
Power (W)	95	
Maximum Vacuum (mbar)	13	
Pump Speed	0 mmHg	25L/min
	100 mmHg	22L/min
	200 mmHg	20L/min
	300 mmHg	12L/min
	400 mmHg	7.5L/min
	550 mmHg	5L/min
Suction Nozzle Size (mm)	5L/min	
Motor Speed (rpm)	1450	
Weight (kg)	7	
Noise (dB)	50	

### Recirculating Chiller Catalogue No: 18200746

Model	ST-F250	
Temperature Range (°C)	-10 ~ +40	
Temperature Stability (°C)	± 0.5	
Temperature display	LED	
Pump capacity flow pressure (bar)	25L/min	
Cooling capacity	100 mmHg	22L/min
	200 mmHg	20L/min
	300 mmHg	12L/min
	400 mmHg	7.5L/min
	550 mmHg	5L/min
Suction Nozzle Size (mm)	5L/min	
Motor Speed (rpm)	1450	
Weight (kg)	7	
Noise (dB)	50	

