

## Water Still

- ST-WS4L Water Still is designed for reliable and continuous operation.
- Chromium plated heater to produce 4L/Hour, Distilled Water, for bench top or wall mounted.



### Specification of Water Still

Model	ST-WS4L
Output	4L/ HR
Conductivity (Raw Water)	3.0 – 4.0 us/cm
Conductivity (Using Treated Water)	1.5 – 2.0 us/cm
pH	5.5 – 6.0
Distillate Quality	Pyrogen Free
Water Supply	Min. flow rate, 1L/min
Heater	3Kw Chromium Plated with thermostat to cut off heater in case of water supply disruption
Glassware	Borosilicate material
Connection	Screw thread connector for condenser and boiler
Dimension (W x D x H)	50W x 15D x 45H cm
Weight	Approx. 12Kg
Voltage	220 – 240V 50Hz, 3Kw