

All Glass Holder



- This product is made of high quality extra hard glass.
- It is sparkling and crystal clear.
- Its thickness is even. The good performance of withstanding voltage and its leakproofness
- The size complies with the international standard.
- It can be used for sterilization at high temperature and high pressure.
- It is mainly used in the filtration of HPLC mobile phase and has the function of outgas to guarantee the clean of mobile phase and avoid the blockage of HPLC liquid flow path.

Specification of All Glass Holder

Part Number	Model	Funnel	Support Screen	Receiving Flask	Suitable Membrane
ST-0202	G1	300 ml	Borosilicate Glass	1000 ml	47mm or 50mm
ST-0202-1	G1-1	500 ml	Borosilicate Glass	2000 ml	58mm or 60mm
ST-0203	G2	300 ml	Borosilicate Glass	1000/500 ml	47mm or 50mm
ST-0207	F300	300 ml	Borosilicate Glass	-	47mm or 50mm

Eyewash Emergency Shower Station

Model: ST-WG7015F

This Floor Mounted Eyewash Emergency Shower combination of self draining shower head and eyewash with foot pedal.

- **Specification for Shower**
Shower Head : Brass
Stand Pipe : Stainless Steel (304 Grade)
Activation : Stainless Steel pull rod operators a stay open stainless steel ball valves
- **Specification for Eyewash**
Bowl : Stainless Steel (304 Grade)
Eyewash Head : Rubber Cups
Activation : Push flag operators a stay open stainless steel ball valves

