



Burners, alcohol, butane, propane and natural gas



ALCOHOL BURNER
Screwcap jar with wick.
Capacity 120 ml.
Part No. **1001256**



ALCOHOL BURNER
Chrome plated top with
wick height control.
With lid and chain.
Capacity 85 ml.
Part No. **1001434**



BUNSEN BURNER
With both gas and air regulator.
Tube Ø exterior: 11 mm. Total height: 140 mm.
Part No.
7000134 For butane/propane gas.
7009133 For natural gas.



SIMPLE BUNSEN BURNER
Air regulator only, nickel plated.
External tube Ø : 11 mm. Total height: 130 mm.
Part No.
7000044 For butane/propane gas.
7000937 For natural gas.



TECLU BURNER
Conforms to CE 370.
With air and gas regulator.
External tube Ø: 13 mm. Total height: 180 mm.
Part No.
7012019 For butane/propane gas.
7094519 For natural gas.



MEKER BURNERS
Conforms to CE 370.
With gas and air regulator.

Model Small. External tube Ø at its maximum point:
22 mm. Total height: 175 mm.
Part No.
7001108 For butane/propane gas.
7001109 For natural gas.



MEKER BURNER
With permanent flame. Gas and air control.
External tube Ø at its widest point: 22 mm.
Total height: 175 mm.
Part No.
7001106 For butane/propane gas.
7001107 For natural gas.



NATURAL GAS BURNERS
Conforms to CE 370.
With air and gas regulator.
With approved gas pressure connection 1/2" nut and
olive, supplied as an accessory.
External tube Ø: 11 mm. Total height: 165 mm.
Part No. For gas Working Normal Weight
pressure power Kg
7001551 natural 18 mbar 0.69 kw 1.6

Model Medium. External tube Ø at its maximum point:
27 mm. Total height: 185 mm.
Part No.
7001110 For butane/propane gas.
7001111 For natural gas.

Model Large. External tube Ø at its widest part:
32 mm. Total height: 190 mm.
Part No.
7001112 For butane/propane gas.
7001113 For natural gas.

ACCESSORY
Gas tube with connections (see page 191).