



Storage refrigerators “Stocklow” L and G

WITH FAN ASSISTED CIRCULATION.

MODEL L, WITH METAL DOOR. TEMPERATURE CONTROLLABLE FROM +1 °C TO +15 °C.

MODEL G, WITH DOUBLE GLAZED ANTI MIST DOOR. TEMPERATURE CONTROLLABLE FROM -2°C TO +16 °C.



APPLICATIONS

Preservation and storage of products that require below ambient temperatures.

FEATURES

- Chemically resistant.
- Tropical hermetic compressor.
- Automatic defrost with electrical water evaporation system.
- Support rack and shelves.
- Adjustable front level supports.
- The door can be hung to open left or right.
- Door lock.

There are two versions:

Stocklow L which is completely made of metal, both interior and exterior, Aisi 310 stainless steel door.

Stocklow G inside made of Aisi 310 stainless steel and epoxy white external case with glazed door that enables internal view. Wheel base support for easy removal.

CONTROL PANEL

- Mains power, ON/OFF switch.
- Electronic temperature regulator.
- Digital temperature display.
- Interior light only (Double glazed door model).



Model Stocklow L, with metallic door.



Model Stocklow G, with metallic door.

MODELS	Part No.	Capacity litres	Door	Height / Width / Depth (interior) cm	Height / Width / Depth (exterior) cm	Number of shelves	Power W	Weight Kg
STOCKLOW L	2101276	600	Metallic	155 51 65	215 70 83	4	250	151
STOCKLOW G	2101277	600	Glass	155 51 65	215 70 83	4	250	151

Comes complete with steel shelves, PVC laminated.



Chest Freezers “Combatemp”

CONTROLLABLE TEMPERATURE RANGE FROM -10 °C TO -28 °C.

CONTROLLABLE TEMPERATURE RANGE FROM -10 °C TO -40 °C.

APPLICATIONS

Suitable for preservation and freezing of sample.

FEATURES

- Epoxy covered external case. Interior of aluminium, with guide supports for baskets or trays.
- The interior of the door is made of easy to clean, strong, chemically resistant, moulded ABS.
- Drop down lockable door with compensating support hinges. Hermetically sealed compressor, silent and vibration free.

- Automatic interior light when the door is opened.
- Digital electronic temperature control.
- N.T.C. Sensor.

CONTROL PANEL

- Temperature control.
- Digital live and preset temperature display.
- Push-button temperature selector.



MODELS	Part No.	Capacity litres	Freezing capacity in 24hrs Kg	Temperature up to °C	Height / Width / Depth (interior) cm	Height / Width / Depth (exterior) cm	Power W	Weight Kg
COMBATEMP S-30	2101040	203	28	-28	55 71 50	86 83 73	160	65
COMBATEMP L-30	2101041	352	38	-28	70 100 52	91 113 76	200	75
COMBATEMP S-40	2101181	203	30	-40	55 71 50	86 83 73	230	70
COMBATEMP L-40	2101182	352	40	-40	70 100 52	91 113 76	400	80

Note: The codes 2101040 and 2101181 are supplied with 3 baskets. The codes 2101041 and 2101182 are supplied with 2 baskets.