

CRUCIBLE TONGS

Original "J.P.Selecta" design with opposite equilibrium.
Made of Electro-bright finished AISI 304 stainless steel.
Part No. **1001330** Total length 400 mm.



CRUCIBLE TONGS

With bow, curved tips and flat hinge.
Electro-bright AISI 304 stainless steel finish.
Part No. **1010420** Total length 200 mm.
Part No. **1010430** Total length 300 mm.
Part No. **1010440** Total length 400 mm.
Part No. **1010450** Total length 500 mm.




CRUCIBLE TONGS

With bow, curved tips.
With thermally protected plastic coated handles.
Part No. **1001590** Total length 220 mm.
Part No. **1001591** Total length 330 mm.




STERILIZATION TONGS

Made of AISI 304 electro-bright finished stainless steel.
Total length 200 mm.
Part No. **1002801** Straight end.
Part No. **1002802** Curved end.
Used for bandages, syringes, instruments, test tubes etc.



TONGS FOR FLASKS AND TUBES

Electro-bright finished AISI 304 stainless steel.
For diameters from 15 to 60 mm.
Tong tips covered with stabilized glass fibre.
Part No. **1028130** Total length 300 mm.
Part No. **1028140** Total length 400 mm.




TONGS FOR FLASKS AND TEST TUBES

Sterilization Tongs. Made of AISI 304 electro-bright finished stainless steel. For diameters from 25 to 60 mm.
Tong ends covered with stabilized glass fibre.
Part No. **1028230** Total length 300 mm.
Part No. **1028240** Total length 400 mm.



DISH TONGS

Made of electro-bright finished AISI 304 stainless steel.
Part No. **1000283** Total length 250 mm.




BEAKER TONGS

Made of AISI 304 electro-bright finished stainless steel.
Made for beakers from 150 to 1000 ml.
For diameters from 60 to 100 mm.
Tong tips covered with stabilised glass fibre.
Part No. **1000284** Total length 450 mm.



BEAKER TONGS "PROTECVAS"

Clamps beakers up to 2000 ml.
With a span from 50 to 120 mm Ø.
Jaws protected with rubber pads.
Total length 300 mm.
Part No. **7000532**



BEAKER TONGS

Original "J.P.Selecta" design with opposite equilibrium.
Made of Electro-bright finished AISI 304 stainless steel.
For beakers from 150 to 1000 ml.
For diameters from 60 to 100 mm.
Tong ends covered with a glass fibre sleeve.
Part No. **1001331** Total length 450 mm.