



LABORATORY SCISSORS

AISI 410 stainless steel.

- 1 Scissors universal "Multi purpose"** with a serrated/ dented blade. Suitable for cutting hard plastic, cardboard and soft wire etc.
Part No. **1001468** Cut length 50 mm. Total length 180 mm.
- 2 Long laboratory scissors.**
Part No. **1001469** Cut length 100 mm. Total length 250 mm.
- 3 Straight multi purpose scissors.** General use and bandage dressing.
Part No. **1001470** Cut length 50 mm. Total length 140 mm.
- 4 Curved multi purpose scissors.** General use and bandage dressing.
Part No. **1000806** Cut length 40 mm. Total length 145 mm.
- 5 Dissection scissors**
Part No. **1000807** Cut length 35 mm. Total length 125 mm.
- 6 Small dissection scissors**
Part No. **1001471** Cut length 25 mm. Total length 110 mm.
- 7 Flat serrated scissors** ideal for laminates, negatives, chromatography, etc.
Part No. **1001472** Flat Length 40 mm. Total length 160 mm.



HOFFMAN CLAMPS / MOHR CLAMPS

- 1 Hoffman clip dual, open-close.**
Part No. **1000063** Maximum opening: 15 mm.
- 2 Hoffman clip closed**
Part No. **1000662** Maximum opening: 15 mm.
- 3 Hoffman clip open**
Part No. **1000661** Maximum opening: 15 mm.
- 4 Hoffman clip dual medium, open-closed.**
Part No. **1001467** Maximum opening: 25 mm.
- 5 Hoffman clip large, open and close.**
Part No. **1000064** Maximum opening: 30 mm.
- 6 Mohr clamp, stainless steel AISI 304 electro bright finish.**
Part No. **1000067** Total length: 70 mm.
- 7 Mohr clamp, nickel wire.**
Part No. **1001337** Total length: 70 mm.



FORCEPS FOR LIFTING SLIDES

Made of AISI 304 stainless steel, with flat jaws.
Length: 105 mm. Suitable for lifting, slides, watch glass or small flat object.
Part No. **1000105**



DISSECTING FORCEPS

- 1 Dissecting Straight forceps round tips, made of AISI304 stainless steel.**
Part No. **1000580** Total length 140 mm.
- 2 Dissecting Straight forceps fine tips, made of AISI304 stainless steel.**
Part No. **1000801** Total length 140 mm.
- 3 Dissecting curved forceps round tips, made of AISI 304 stainless steel.**
Part No. **1000865** Total length 140 mm.
- 4 Dissecting curved forceps fine tips, made of AISI304 stainless steel.**
Part No. **1000864** Total length 140 mm.

MAXIMUM OPENING

- 5 Sprung straight forceps with fine tips, made of AISI 304 stainless steel.**
Part No. **1001189** Total length 160 mm.
- 6 Sprung straight forceps with curved fine tips, made of AISI 304 stainless steel.**
Part No. **1001190** Total length 160 mm.
- 7 Curved point, AISI 304 stainless steel forceps.**
Part No. **6000128** Total length 160 mm.
- 8 Fine point, AISI 304 stainless steel forceps.**
Part No. **1001191** Total length 120 mm.



HEAT AND SPILLS PROTECTION PAD GLOVE

Made of flexible silicon that fits all hand sizes.
Part No. **1001255**