

SUPPORT STAND KITS



700 high x 700 mm long.

Made up of:
8 rods 700 mm long.
2 rods 350 mm long.
16 "P" type bossheads.
4 foot supports.
1 Allen key.

Part No. **600090** Dur-aluminium

Part No. **600190** AISI 304 stainless steel




1000 high x 1000 mm long.

Made up of:
10 rods of 1000mm long.
2 rods 700 mm long.
3 rods 350 mm long.
34 "P" type bossheads.
6 foot supports.

1 Allen key.

Part No. **6000902** Dur-aluminium

Part No. **6001902** AISI 304 stainless steel




700 high x 1000 mm long.

Made up of:
4 rods 1000 mm long.
5 rods 750 mm long.
2 rods 350 mm long.
20 "P" type bossheads.
4 foot supports.
1 Allen key.

Part No. **6000901** Dur-aluminium

Part No. **6001901** AISI 304 stainless steel



1000 high x 2000 mm long.

Made up of:
5 rods 2000 mm long.
7 rods 1000 mm long.
8 rods 200 mm long.
2 rods 350 mm long.
45 "P" type bossheads.
8 foot supports.

1 Allen key.

Part No. **6000903** Dur-aluminium

Part No. **6001903** AISI 304 stainless steel



SUPPORT "STAND"

Solid base, three vertical rods of 800 mm in height, two horizontal rods of 500 mm in length, six "P" type bossheads. Made of Dur-aluminium.

Weight: 6.5 Kg.

Part No. **6000482**

CENTRAL ROD RETORT STAND

Vertical rod of Dur-aluminium of 10 Ø x 600mm high.

Weight: 3 Kg.

Base dimensions 315 x 200 mm.

Part No. **6000046**



TRIPOD SUPPORT "PASS"


Original Phywe.

To support square and round rods from 6 to 14 mm dia.

Weight: 2 Kg.

Base with levelling screws.

Part No. **6200255**



DOUBLE BASE SUPPORT T "MINI"

Small foot support with screw clamp for removable rods from 8 to 12 mm Ø.

Weight: 1.2 Kg.

Part No. **6001103**



STANDARD RODS

Made from AISI 340 stainless steel. Suitable for double "T MINI" and Tripod "PASS"

Part No. **6001160** 8 Ø x 500 mm length

Part No. **6000534** 10 Ø x 600 mm length

Part No. **6000535** 12 Ø x 700 mm length