

## Digital Ultrasonic Cleaner With Heating

- Double wall construction, made entirely of stainless steel.
- LED display of timer, up to 80 minutes and temperature up to 90°C.
- Dual-band frequency, adjustable, 40/59 kHz with or without drainage built-in.
- Supply complete with basket and lid.



### Specification of Digital Ultrasonic Cleaner With Heating

Model	ST-UB2200LDT	ST-UB3200LDT	ST-UB3300LDT	ST-UB5200LDT	ST-UB7200LDT	ST-UB8200LDT
Capacity	3L	4L	6L	10L	22L	22L
Frequency	40/59 kHz					
Power	120W	150W	180W	250W	500W	600W
Drainage	No	No	Yes	Yes	Yes	Yes
Timer	Yes	Yes	Yes	Yes	Yes	Yes
External Dimension (mm)	260 x 160 x 290	320 x 175 x 290	350 x 175 x 310	350 x 264 x 310	550 x 320 x 310	550 x 320 x 310
Internal Dimension (mm)	240 x 140 x 100	300 x 150 x 100	300 x 150 x 150	300 x 240 x 150	500 x 300 x 150	500 x 300 x 150

## Ultrasonic Cleaner

- Double wall construction, made entirely of stainless steel.
- Draining stopcock only for model 6 litres onwards.
- Timer 0-120 minutes.
- 40kHz sound wave frequency for model 2L, 3L, 4L, 6L, 10L & 22L, 40 or 40/59kHz & 28/40/59 kHz(factory preset) for larger model.
- Supply complete with basket and lid.



### Specification of Digital Ultrasonic Cleaner With Heating

Model	ST-UB80D	ST-UB120D	ST-UB150D	ST-UB180D	ST-UB250D	ST-UB300D	ST-UB500D	ST-UB600D
Capacity	2L	3L	4L	6L	10L	10L	22L	22L
Frequency	40 kHz				40/59 kHz		28/40/59 kHz	
Power	80W	120W	150W	180W	250W	300W	500W	600W
Drainage	No	No	No	Yes	Yes		Yes	
Timer	Yes	Yes	Yes	Yes	No		No	
External Dimension (mm)	173 x 160 x 190	260 x 160 x 240	320 x 175 x 240	320 x 175 x 310	350 x 265 x 310	350 x 265 x 310	550 x 325 x 310	550 x 325 x 310
Internal Dimension (mm)	150 x 140 x 100	240 x 140 x 100	300 x 150 x 100	300 x 240 x 150	300 x 240 x 150	300 x 240 x 150	500 x 300 x 150	500 x 300 x 150