

Ductless Fume Hood

- SASTEC Ductless fume hood provide protection to the operator and the environment.
- Chemical fume and gasses are being filtered and contained in the carbon filter.
- Clean air are being re-directed back into the laboratory.
- The inflow air with an average 0.5 m/s passes through the front to the work area.
- Negative pressure is maintained to ensure no chemical fume or gasses escape from the work area.
- Standard supply comes with Organic Solvent Carbon Filter.
- Optional filter is available for Acid Fumes and Safety Ulpa Filter. (organic solvent + acidic fume + ulpa filter)
- The hood is constructed with EG steel (Electro Galvanized) for 1mm thickness come with epoxy coating and work top is made by epoxy resin.
- Model 4FT, 5FT and 6FT comes with base support.
- MoModel 3Ft, 4Ft, 5Ft & 6Ft all surround view (transparent).
- Optional: UV Lamp
- Accessories:
 - * LED Lamp
 - * Water Tap with Oval Sink
 - * Gas Tap
 - * Base Cabinet
 - * Internal Phenolic Resin Panel
 - * Vaneometer for airflow indication



Specification of Ductless Fume Hood

Model	ST-DF3	ST-DF4	ST-DF5	ST-DF6
Voltage	200~240V 50Hz			
Power (W)	150W	150W	150W x 2	150W x 2
Noise Level	< 58dB(A)	≥58 dBa	< 60dB(A)	< 60dB(A)
Control System	Rocket Switch			
Fluorescent Light	Intensity >500 Lux			
Working Zone (W x D x H) mm	822 x 744 x 585	1180 x 778 x 745	1480 x 778 x 745	1780 x 778 x 745
Overall Dimension (W x D xH) mm	896 x 779 x 1035	1280 x 818 x 2300	1580 x 818 x 2300	1880 x 818 x 2300
Weight	45 kg	75 kg	95 kg	115 kg
Construction	Top Hood	Chemical Resistant Paint & E.G Steel		
	Hood Wall	Polycarbonate Clear Plane		
	Work Base	Epoxy Resin		
Air Flow System	Inflow Velocity	0.5 m/s ±0.05		
	Indiside Volume of Hood	0.383 m ³	0.571 m ³	0.716 m ³
	Volume of Treated Air (hour)	930 m ³ /hour	1270 m ³ /hour	1609 m ³ /hour
3 Type Filter in 1 Hood	Volatile Organic Compound (VOC)	Yes		
	Acidic Gas (Acid)	Optional		
	ULPA Filter	Optional		
Stand	Frame c/w locking castle wheel	Optional	Yes	Yes
	Overall Dimension c/w frame	896 x 779 x 1741	N/A	N/A
	Weight	55 kg	N/A	N/A
SASTEC Touch Pad	Sash	Motorized Up & Down		
	Power	On/Off		
	Lamp	On/Off		
	Blower	On/Off		
Cabinet	Cabinet	Optional		
	Weight	65kg	105kg	125kg
	Overall Dimension c/w Frame (W x D xH) mm	896 x 779 x 1035	1280 x 818 x 1906	1580 x 818 x 1906

Diagram & Part List

- ① TUBULAR MOTOR
- ② BLOWER
- ③ ROCKET SWITCH
- ④ FLUORESCENT LAMP
- ⑤ SOCKET OUTLET SINGLE
- ⑥ POWER CORD
- ⑦ FRAME C/W CASTER WHEEL

