

Top Loading Autoclave (Vertical Type)

- SASTEC Microprocessor Controller (Standard Model)
 - Fully automatic operation
 - LCD Screen Display of working parameters.
 - Analog pressure gauge
 - Safety Valve for over pressure, for low water level.
 - Safety interlocking prevent opening of cover during sterilization cycle.
 - Alarm for end of cycle, safe to open over.
 - Single one action door lock and release.
 - Stainless Steel 316 heating element to resist corrosion
 - Stainless Steel 304 internal chamber and external body.
 - Conforming to internation pressure regulation and safety requirement.
- Optional features
- : Vaccum Drying Function
 - : Printing
 - : Additional Water Reservoir
- With JKPP Certification.
 - With TUV Certification.



Specification of Top Loading Autoclave (Vertical Type)

Model	ST-TLA40	ST-TLA60	ST-TLA80
Capacity (L)	53L	70L	90L
Chamber Dimension (mm)	370 x 380 mm	370 x 560	370 x 750
External Dimension (mm)	710W x 540L x 880H mm	710W x 540L x 980H mm	710W x 540L x 1113H mm
Gross Weight (Kg)	95 Kg	105 Kg	120 Kg
Voltage (Volt)	220 Volt		
Power (KW) 15amp	3.5		
Temperature °C	50 ~ 134°C		
Timer	0 -500 Minutes		
Accessories included Platform	1		
Wire Baskets	1	2	3
With Vacuum	ST-TLA40V	ST-TLA60V	ST-TLA80V
With Inner Reservoir	ST-TLA40S	ST-TLA60S	ST-TLA80S