

Upright Two Door Freezer



- Temperature Range- Freezer : +2 to +8°C/ -10 to -40 °C
- Microprocessor temperature controller, LCD display
- Audio/Visual alarm for High/Low temperature, power failure, door open and sensor failure.
- 90mm CFC-free insulation
- Danfoss compressor and EBM motor fan.
- Double sealing with lock
- Painted metal external and ABS material internal.
- Three shelves for upper and lower chamber.
- Mounted on castors.
- CFC-free refrigerant R134a .
- Certification: ISO 9001, ISO 14001, ISO 13485, TUV
- Power supply: 220V 50Hz (110V 60Hz is available)
- Warranty period: 1 year on all parts.

Specification of Up Right Fridge-Freezer

Model	ST-FF40-253	
Type	Upper	Lower
Capacity (L)	154L	99L
External Dimension (W x D x H) mm	620 x 670 x 1780	
Electrical Power	220V,50/ 60HZ, 110V,50/60HZ	
Input Power (W)	400W	
Refrigeration System	International famous compressor and Germany EBM fan motor	
Controller	Microprocessor Control with LCD Displayer	
Alarms	High&low temperature alarm, Power failure alarm, Door open alarm, sensor failure alarm	
Optional	Remote alarm system	
Weight	74Kg	